

**CONFIDENTIAL**

**THE UNITED REPUBLIC OF TANZANIA  
NATIONAL EXAMINATIONS COUNCIL OF TANZANIA  
CERTIFICATE OF SECONDARY EDUCATION EXAMINATION 2023**

**032/2A**

**CHEMISTRY 2A  
ACTUAL PRACTICAL A**

**3 HOURS PRACTICAL ADVANCE INSTRUCTIONS**

**1.0 IMPORTANT**

- 1.1 GREAT CARE MUST BE TAKEN NOT TO DIVULGE THESE INSTRUCTIONS TO BOTH CANDIDATES AND UNAUTHORIZED PERSONS EITHER DIRECTLY OR INDIRECTLY.
- 1.2 MAKE SURE THAT THE CANDIDATES ARE PROVIDED WITH CHEMICALS AND APPARATUSES AS INDICATED IN THESE PRACTICAL INSTRUCTIONS ONLY AND NOT OTHERWISE.

**2.0 PREPARATION AND LABELLING OF CHEMICALS AND APPARATUSES**

**2.1 Question 1**

- Prepare 0.1 M HCl, label it V1 and allow 125 cm<sup>3</sup> per candidate. Choose from the following table the appropriate volume of stock (concentrated) solution to be diluted to make 1 L solution basing on the specification indicated on the bottle of stock solution.

S/n	Percentage Purity/Assay or its Average when Given in a Range	Density of a Stock Solution (g/cc)	Volume of Stock Solution to be Diluted to make 1 L (cm <sup>3</sup> )
1	31 - 32	1.16	9.81
2	≈ 34	1.18	9.09
3	≈ 35	1.18	8.83
4	≈ 36	1.18	8.58
5	≈ 37	1.18	8.35
6	≈ 38	1.18	8.13

Dilute the selected volume with water to make 1 L of 0.1 HCl solution.

- Prepare 0.1 M of sodium hydroxide by dissolving 4 g of the base in water to make 1 L solution. Label it V2 and allow 100 cm<sup>3</sup> per candidate.
- Provide each candidate with methyl orange (MO) indicator.
- Provide each candidate with two test tubes.
- Provide blue and red litmus papers, each one per candidate.

## CONFIDENTIAL

- Provide each candidate with 1 burette, 1 titration flask (250 cm<sup>3</sup>), 1 pipette (20 cm<sup>3</sup> or 25 cm<sup>3</sup>), 1 white tile, 1 glass/plastic funnel, 1 pipette filler, 2 test tubes, 1 retort stand and a clamp.

### 2.2 Question 2

- Prepare 0.2 M Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>·5H<sub>2</sub>O by dissolving 49.636 g of hydrated sodium thiosulphate in water to make 1 L solution. Label it RR and allow 70 cm<sup>3</sup> per candidate.
- Prepare 1 M HCl, label it SS and allow 70 cm<sup>3</sup> per candidate. Choose from the following table the appropriate volume of stock (concentrated) solution to be diluted to make 1 L solution basing on the specification indicated on the bottle of stock solution.

S/n	Percentage Purity/Assay or its Average when Given in a Range	Density of a Stock Solution (g/cc)	Volume of Stock Solution to be Diluted to make 1 L (cm <sup>3</sup> )
1	31 - 32	1.16	99.78
2	≈ 34	1.18	90.88
3	≈ 35	1.18	88.28
4	≈ 36	1.18	85.83
5	≈ 37	1.18	83.51
6	≈ 38	1.18	81.31

Dilute the selected volume with water to make 1 L of 0.1 M HCl solution.

- Provide each candidate with a stop watch and a thermometer.
- Provide each candidate with 60 cm<sup>3</sup> distilled water.
- Using a blue pen, draw a uniform mark M on one side of a 8 cm × 8 cm piece of white paper and allow one piece per candidate.
- Provide each candidate with a 50 cm<sup>3</sup> beaker.
- Provide each candidate with a 10 cm<sup>3</sup> measuring cylinder.